EPA LIST N

EPA approved on List N, for use against SARS-CoV-2, the cause of COVID-19.



CLEANCIDE GERMICIDAL PRODUCTS

A SMARTER, SAFER WAY TO CLEAN + DISINFECT

Citric acid-based CleanCide germicidal wipes deliver powerful efficacy against pathogens — without the harmful, irritating effects of other chemical disinfectants — in an easy-to-use wipe format that enables one-step, no-rinse cleaning and disinfection.



POWERFUL DISINFECTION - BROAD EFFICACY

Proven, broad spectrum, hospital-grade antimicrobial efficacy — virucidal, germicidal, fungicidal and tuberculocidal. Meets OSHA requirements for bloodborne pathogens and EPA List N criteria for use against SARS-CoV-2, the virus that causes COVID-19.



IDEAL FOR HARD, NON-POROUS SURFACES

Formulated for fast, effective cleaning and disinfection of hard, nonporous surfaces such as stainless steel, chrome, glass, plastic and vinyl. Ideal for use on door knobs, bed rails, monitoring equipment and other digital monitors, telephones, light switches, faucet handles and other high-touch surfaces.



CITRIC ACID-BASED - NON-IRRITATING

A better alternative to quat-, alcohol- or bleach-based disinfectants. EPA Category IV (lowest level of toxicity) requiring no precautionary statements or PPE — will not irritate skin, eyes or lungs. Highly biodegradable ingredients, leaving no harmful or irritating residues.



ONE-STEP, NO-RINSE CLEANING & DISINFECTION

Proven efficacy in the presence of 5% organic soil load. 5-minute contact time for disinfection. Enables one-step cleaning and disinfection — no pre-cleaning or post-rinsing steps required.







PRODUCTS

PRODUCT #	DESCRIPTION	PACKAGING
2130-02	CleanCide Ready-To-Use*	12 Quarts
2130-07	CleanCide Ready-To-Use*	5 Gal Pail
2130-09	CleanCide Ready-To-Use*	55 Gal Drum
3130C-160-CS	CleanCide Wipes – 160 Count Canister	12 Canisters
3130FP-80-CS	80 count flow pack	12 Packs
3130B-400D	CleanCide Wipes – 400 Count Tub	1 Tub



¹Sanitizing claim only applies to RTU liquid







CLEANCIDE EFFICACY DATA

FPA REGISTRATION NUMBER 34810-35

DISINFECTING Bactericidal activity (in the presence of 5% blood serum) 5 MINUTES ATCC Carbapenem-resistant Klebsiella pneumoniae ATCC BAA-1705 Enterococcus faecalis Vancomycin Resistant VRE ATCC 51575 Escherichia coli O157:H7 ATCC 43888 ATCC 19117 Listeria monocytogenes Pseudomonas aeruginosa ATCC 15442 Salmonella enterica (formerly Salmonella Choleraesuis) ATCC 10708 Staphylococcus aureus ATCC 6538 Staphylococcus aureus - MRSA ATCC 33592 Staphylococcus epidermidis - MRSE ATCC 51625 In the presence of 5 % blood serum, is effective against 3 MINUTES STRAIN AND/OR ATCC Hepatitis B (HBV), Duck Hepatitis B (DHBV) as surrogate Grimaud Strain Hepatitis C (HCV), Bovine Viral Diarrhea Virus (BVDV) as surrogate Strain NADL Strain 1059 Rhinovirus **5 MINUTES** ATCC VR-846 Adenovirus Type 2 Feline Calicivirus (surrogate for Norovirus) Strain F-9 Herpes Simplex Virus Type 1 ATCC VR-733 ATCC VR-734 Herpes Simplex Virus Type 2 Human Immunodeficiency Virus Type 1 (HIV-1) (AIDS Virus) Strain HTLV-IIIB Influenza A Virus Strain Hong Kong H3N2 Strain Pandemic 2009 H1N1 Influenza A Virus (formerly called swine flu) Respiratory Syncytial Virus ATCC VR-26 Rotavirus Strain WA, Univ. of Ottawa ATCC VR-119 Vaccinia Virus Fungicidal activity (in the presence of 5% blood serum) ATCC 10231 Candida albicans Trichophyton interdigitale (T. interdigitale) ATCC 9533 Tuberculocidal activity at 20°C/68°F on previously cleaned hard, non-porous, inanimate surfaces **5 MINUTES** STRAIN Mycobacterium bovis (Mycobacterium tuberculosis) (M. bovis) (TB) Strain BCG (Organon Teknika Corp.) SANITIZING* 60 SECONDS ATCC Klebsiella pneumoniae ATCC 4352 FOOD CONTACT SANITIZER: 60 SECONDS ATCC Staphylococcus aureus ATCC 6538 Escherichia coli ATCC 11229



AUTHORIZED DISTRIBUTOR

CLEANER SOLUTIONS INTERNATIONAL, LLC.





^{*} Efficacy claims apply to RTU product only